

acetate ($\text{CH}_3\text{COONa}\cdot 3\text{H}_2\text{O}$) and again stir until dissolved. Dilute to 500 mL with distilled water. This solution is stable for several weeks if stored in the dark.

- 6.3 Nitrite stock solution: 1.0 mL = 0.10 mg $\text{NO}_2\text{-N}$. Dissolve 0.1493 g of dried anhydrous sodium nitrite (24 hours in desiccator) in distilled water and dilute to 1000 mL. Preserve with 2 mL chloroform per liter.
- 6.4 Nitrite standard solution: 1.0 mL = 0.001 mg $\text{NO}_2\text{-N}$. Dilute 10.0 mL of the stock solution (6.3) to 1000 mL.

7.0 Procedure

- 7.1 If the sample has a pH greater than 10 or a total alkalinity in excess of 600 mg/L, adjust to approximately pH 6 with 1:3 HCl.
- 7.2 If necessary, filter the sample through a 0.45 μ pore size filter using the first portion of filtrate to rinse the filter flask.
- 7.3 Place 50 mL of sample, or an aliquot diluted to 50 ml, in a 50 mL Nessler tube; hold until preparation of standards is completed.
- 7.4 At the same time prepare a series of standards in 50 mL Nessler tubes as follows:

mL of Standard Solution 1.0 mL = 0.001 mg $\text{NO}_2\text{-N}$	Conc., When Diluted to 50 ml, mg/L of $\text{NO}_2\text{-N}$
0.0	(Blank)
0.5	0.01
1.0	0.02
1.5	0.03
2.0	0.04
3.0	0.06
4.0	0.08
5.0	0.10
10.0	0.20

- 7.5 Add 2 mL of buffer-color reagent (6.2) to each standard and sample, mix and allow color to develop for at least 15 minutes. The color reaction medium should be between pH 1.5 and 2.0.
- 7.6 Read the color in the spectrophotometer at 540 nm against the blank and plot concentration of $\text{NO}_2\text{-N}$ against absorbance.

8.0 Calculation

- 8.1 Read the concentration of $\text{NO}_2\text{-N}$ directly from the curve.
- 8.2 If less than 50.0 mL of sample is taken, calculate mg/L as follows:

$$\text{NO}_2 - \text{N, mg/L} = \frac{\text{mg/L from std. curve} \times 50}{\text{mL sample used}}$$

9.0 Precision and Accuracy

9.1 Precision and Accuracy data are not available at this time.

Bibliography

1. Standard Methods for the Examination for Water and Wastewater, 14th Edition, p 434, Method 420, (1975).